







Aliya[®] PEF System

For Soft Tissue
Ablation



Aliya[®] PEF System


a fundamentally
different approach to
soft tissue ablation



Cell Homeostasis Disruption

The Aliya PEF system is a revolutionary soft tissue ablation device that delivers high voltage, short duration electrical energy to disrupt cell homeostasis.


The proprietary Aliya PEF mechanism creates a cascade of biochemical processes that results in **regulated cell death**, a nonthermal mechanism of action that offers unique procedural and patient impacts.



**Aliya PEF's
Ablation
Zone**

Intact Airway

Central tissue void
due to punch biopsy



Favorable Safety & Predictability*

- Provides a predictable, consistent ablation zone unaffected by heat sinks
- Uses a small 19G needle to deliver Aliya PEF ablation energy
- Preserves extracellular matrix and surrounding tissue
- Has been used safely near airways, vessels and pleural surfaces

Clinical data suggests a procedural burden similar to diagnostic biopsy



Antigen Release

Potential for Immune Activation*


- Works with the body's biology to achieve durable target tissue ablation
- May promote antigen release that could activate the immune system

Adjunct to Standard of Care*


- Does not interfere with treatment plan and sequencing post-PEF
- Can be performed at the same time as diagnostic biopsy

Does not interfere with standard of care therapies

Aliya Electrode



Aliya Needle



- Predictable ablation zone
- 10 mm - 13 mm

ALIYA® PULSED ELECTRIC FIELD SYSTEM EQUIPMENT

Part Number	Description	Qty	UOM
GTI-00017	Aliya Pulsed Electric Field System Generator 2.5	1	Each
GTI-00009	Cardiac Monitor, Synchronization, IVY Model 7600	1	Each
GTI-00053-00	Aliya 2.5 System 1 Year Extended Warranty	1	Year

ALIYA® PULSED ELECTRIC FIELD SYSTEM CONSUMABLES

Part Number	Description	Qty	UOM
GTI-00050-5PK	Aliya PEF Ablation Pack (Aliya Needle 19G, Aliya Electrodes)	1	Box (5)
GTI-00023-1PK	Aliya PEF Needle, 19 Gauge	1	Box (5)



Galvanize Therapeutics, Inc.
3200 Bridge Parkway
Redwood City, CA 94065

www.galvanizetx.com

info@galvanizetx.com

[in galvanize-therapeutics](https://www.linkedin.com/company/galvanize-therapeutics)
(650) 268-4252

* References on file.

The Aliya® System is 510(k) cleared in the United States for the surgical ablation of soft tissue. It is not currently commercially available in any other jurisdiction. Galvanize Therapeutics does not promote the off-label use of its products and nothing herein is intended to promote an off-label use of the Aliya System. The Aliya System is a tool for the surgical ablation of soft tissues, and is not intended to treat, cure, prevent or mitigate any specific disease or condition.

CAUTION: Federal (US) law restricts this device to sale by or on the order of a physician.

Important information: Prior to use, refer to the Instructions for Use that are supplied with this device for indications, contraindications, side effects, suggested procedure, warnings and precautions. Galvanize Therapeutics and Aliya are trademarks and may be registered in the US and/or in other countries.